


Life Science

It's in your hands.

www.cryotherm.de | www.cryotherminc.com





Prolog

Our Life Science Portfolio.

Cryopreservation using cryogenic liquefied nitrogen is an established and reliable practice in science, research, medicine, environmental research and engineering. The safe cooling down to -190°C in the gas phase as well as the inert atmosphere allow for optimal freezing and storage of valuable medical and biological samples and materials.

Cryotherm, a leading supplier of storage, transfer and transport systems, has been a renowned manufacturer of cryogenic technology products and services for more than 55 years.

BIOSAFE[®] by Cryotherm has long been synonymous with superior quality and safety. Experience, technology and our know-how make it easy for you to decide on BIOSAFE[®]. Advancing into the future with recent cutting-edge developments, we have expanded our BIOSAFE[®] line-up with new life science products.

Our BIOSAFE®-systems consist of:

- BIOSAFE[®] Cryogenic storage vessels; to freeze and store large quantities of samples.
- BIOSAFE[®] smart, a level controller device and container monitoring system for reliable supply of cryogenic liquefied nitrogen
- BIOSAFE® VENUS an aluminum container series for small to medium sample volumes.
- BIOSAFE[®] MARATHON a dry shipping container for safe transport of your biological samples.

www.cryotherm.de | www.cryotherminc.com

Table of contents

APOLLO® to VENUS.

BIOSAFE[®] - Systems An overview of our life science family.

Safety equipment Your safety is our priority.

BIOSAFE® 120 - 220 - 420 The secure BIOSAFE® wide neck storage vessel series. (for up to 23,100 cryogenic vials 2 ml)

BIOSAFE® 500 - 800 - 2400 The narrow neck BIOSAFE® Stainless Steel Series. (for up to 120,000 cryogenic vials 2 ml)

BIOSAFE[®] smart The 'digital brain' of your biobank.

BIOSAFE[®] iTrack Your intelligent sample management.

BIOSAFE® VENUS The safe storage for your samples. (for up to 6,000 cryogenic vials 2 ml)

BIOSAFE® MARATHON For safe shipment of your samples. (for up to 612 cryogenic vials 2 ml)

BIOSAFE® ALPHA Unpressurized vessel series for transport. (10 - 25 l capacity)

APOLLO[®] Vacuum super-insulated storage and transport vessel for liquid (50 - 350 | capacity)



	2 - 3
	6 - 7
	8 - 10
	11 - 13
	14 - 16
	17 - 19
	20 - 24
	25 - 27
	28 - 29
l nitrogen.	30 - 33



Safety equipment

Your safety is our priority.



Cryo-Gloves, waterproof				
Description Hand protection, 5 fingers, down to -160° C	Length	Size	ltem No.	
Length: ca. 300 mm 1	Wrist	S	78400613	
Length: ca. 300 mm	Wrist	М	78213086	
Length: ca. 300 mm	Wrist	L	77031450	
Length: ca. 300 mm	Wrist	XL	78400823	
Length: ca. 400 mm 💈	Forearm	S	78400843	
Length: ca. 400 mm	Forearm	Μ	78400180	
Length: ca. 400 mm	Forearm	L	78400887	
Length: ca. 400 mm	Forearm	XL	78400856	

Elbow

Elbow

Elbow

M 78400899

L XL

77031451

78211327





Body Protecti	on		
Cryo-Aprons		Beschreibung	Artikelnr.
length: approx. 920 mm, for and splashes for temperatu	re range to -160 °C		77031519
length: approx. 1380 mm, for and splashes for temperatu	77031520		
Cryo-Protective Shield an	d Goggles		
cryo-protective shield and inspection Glass	head protection sh zed fibre foam padde (WxH)		77031519
Schild			
"Wear eye protection"	Adhesive Foil Ø	20 cm Symbol	78210204
"Wear hand protection" 7	Adhesive Foil Ø	20 cm Symbol	78210205







4

Length: ca. 500 mm

Length: ca. 500 mm

Length: ca. 500 mm









BIOSAFE® 120 - 220 - 420

Our secure BIOSAFE® wide neck storage vessel series.

All Advantages at a Glance



Our BIOSAFE® systems 120 to 420 are optimized for maximum capacity with minimized external dimensions. The large, thermodynamically optimized, vacuum-insulated lid that comes with an attractive design makes it possible to access stored materials in just seconds.

Filling of samples takes place through the wide neck opening in the freely accessible storage compartment. The solid construction, made of polished, stainless steel, ensures maximum life time and a low evaporation rate. A combined vacuum closure and safety device for the vacuum chamber provides additional security.

BIOSAFE® 120 - 220 - 420

Our secure BIOSAFE® wide neck storage vessel series.



* available for BIOSAFE[®] 220 / 420





Double the safety!

Mechanically and (optional) electronically lockable lid

Well-thought-out:

Depending on the vessel size, high nitrogen reserves below the shelf floor will last several days.

High storage capacity

for up to 23,100 samples, 2 ml vials

Smooth surfaces, that are easy to disinfect - even on top of the lid.



BIOSAFE® 120 - 220 - 420

Ideal for routine work with biological-medical samples.

Storage Capacities

	~			
Wide-neck BIOSAFE®		120	220	420
Geometric content	[]	151	251	420
2 ml vials	[max./up to]	7.000	13.000	23.100
Baxter 4R9951 bag 50 ml	[max./up to]	288	512	832
and 2 ml vials				1.500
Baxter 4R9955 bag 500 ml	[max./up to]	144	240	480
and 2 ml vials		500	1.500	500
Hemofreeze Gambro Z2003 bag 500 ml	[max./up to]	75	150	270
and 2 ml vials		750	500	500

Technical Data		1		
Wide-neck BIOSAFE®		120	220	420
Neck diameter	[mm]	514	664	854
LIN Capacity under platform	[]	21	35	57
Weight, empty	[kg]	100	175	225
Diameter, outside	[mm]	560	710	900
Diameter, inside	[mm]	514	664	854
Total height, closed	[mm]	1.045	1.070	1.070
Total height, opened	[mm]	1.140	1.180	1.190
Installation altitude	[mm]	625	625	625
Caster diameter	[mm]		80	80
Evaporation rate	[%/day]	2,3	2,0	1,5









BIOSAFE® 500 - 800 - 2400

Our Narrow Neck BIOSAFE[®] series



Our BIOSAFE® systems 500 to 2400 are optimized for maximum capacity with minimum consumption. In the range of 500 liters to 2400 liters and more, the containers can be easily loaded via the side-mounted vessel neck and the smoothrunning carousel. Smooth-running castors ensure the mobility of our BIOSAFE[®] storage vessels and provide even more flexibility.








BIOSAFE® 500 - 800 - 2400

Our Narrow Neck BIOSAFE® Stainless Steel Series.

All Advantages at a Glance



BIOSAFE[®] 500 - 800 - 2400

Our Narrow Neck BIOSAFE® Stainless Steel Series.

Storage capaciti	es	$\overline{\forall}$	7			
Narrow neck BIOSAFE®		-		500	800	2400
Geometric content		[]		490	960	2.400
2 ml vials	[max	./up to]	28	.280	43.120	120.000
Baxter 4R9951 bag 50 ml	[max./up to]		1	.080	on request	on request
and 2 ml vials				650	on request	on request
Baxter 4R9955 bag 500 ml	[max./up to]			480	on request	on request
and 2 ml vials			1	.300	on request	on request
Hemofreeze Gambro Z2003 bag 500 ml	[max	./up to]		288	on request	on request
and 2 ml Ampullen			1	.300	on request	on request

	Technical Data		1		
	Narrow neck BIOSAFE®		500	800	2400
	Neck diameter	[mm]	450	520	830
	LIN Capacity under platform	[]	50	117	300
	Weight, empty	[kg]	320	570	1.300
	Diameter, outside	[mm]	890	1.058	1.650
	Diameter, inside / carousel	[mm]	795	966	1.500
	Total height, closed	[mm]	1.014	1.380	2.100
	Total height, opened	[mm]	1.314	1.800	2.100
	Installation altitude	[mm]	800	800	915
	Caster diameter	[mm]	80	80	100
	Evaporation rate	[%/day]	1,2	0,9	0,6
1	BOGAFE"			BOSHT	-
a.	11				





13

BIOSAFE® smart

BIOSAFE® smart

Smart. Stable. Steady.

Biobanking 4.0

-190°C

660

Big Data

From monitoring to recording to evaluation of the sample cycle - Let data speak.

Touch Screen

Intuitive multifunctional 8.8" panel with colored touch display. Bright and functional.

Conversion Kit*

Your existing BIOSAFE® repository freezer can be updated with new BIOSAFE® smart. Item no. 78221956

* Cryotherm verifies feasibility.

Best durability

Protected by a sturdy glass panel, ideal for aseptic and medical environments.

Select your Software-Package

Smart Config-Version On-site BIOSAFE® monitoring via color display (included).

Smart + Info.View On-site and remote monitoring of BIOSAFE® including secure data analysis and transfer, with the ability to connect up to 254 BIOSAFE® units in up to 10 cryobanks. optional: Item no. 78216677

Connect your Biobank, be it in Berlin, Atlanta or Shanghai...

www.cryotherm.de | www.cryotherminc.com



100 % in control

Thanks to optional electromechanical locking, your samples are protected against access. Twofold access control via RFID transponder and PIN code. Logs all activities according to GMP FDA CRF 21 part 11.



Extensive special equipment

Various options for linking components such as gas bypass, UPS battery, lid lock, additional temperature sensors and much more.

> Smart + messenger Info.Alarm Includes all features of Smart + including Messenger function via encrypted SMS and / or eMail.

optional: Item no. 78216678




BIOSAFE® smart

Smart. Stable. Steady.

The one-stop Biobank





BIOSAFE® iTrack

Your intelligent sample management.

Tracking your Cryogenic Samples

Storage

With BIOSAFE® iTrack you give your samples a digital CV, tracking every step from initial storage to removal of your biological samples in your BIOSAFE®.

iTrack PREMIUM The good-is-better package

iTrack samples DataBase Explorer as with Web Access (interconnectivity) and Sample Data Management

- for integration of 1D 2D scanners 2D flatbed scanners
- for up to five users

Item number: 78216918

- Support during the contract period during normal business hours
- support
- reports

One-time installation fee

ilnstallation fee	Ins
ltem number: 78216916	Iter
Annual service charge (optional)	An


nual service charge (optional) Item number: 78216918



🔸 🗕 🖣 Made in Germany









(f) (in) 🕑 💘 💌

BIOSAFE® smart

Smart. Stable. Steady.

Functional favorite



BIOSAFE[®] smart is an on-site vessel monitoring system with the following abilities:

- Level control and monitoring, automatic filling and temperature display
- Monitoring and evaluation of alarms and events
- Data analysis and transfer according to GMP FDA CRF 21 Part 11
- Double access control via RFID transponder and PIN code
- Temperature monitoring with up to 3 temperature sensors
- Visual and audible alarm on the display
- Forwarding of alarms and temperature
- Other optional equipment: battery backup, warm gas by-pass, software with alarms and messenger

BIOSAFE® smart

X	Item Number	Kurztext St	k
\Box	78221903	BIOSAFE [®] smart120	_
\Box	78221906	BIOSAFE [®] smart220	
\Box	78221909	BIOSAFE [®] smart420	
\Box	78221912	BIOSAFE [®] smart500	
\Box	78221918	BIOSAFE [®] smart1000	_
	78221924	BIOSAFE [®] smart2400	


Software Options					
78216677	BIOSAFE [®] smart+ (Info.View)				
	BIOSAFE [®] smart+ messenger				
78216678	(Info.Alarm)				


O V	Medical Device according to MDD 93/42/EWG						
	Ø	Item Number	Kurztext	Stk			
	\Box	78221879	BIOSAFE®smart120 MD				
	\Box	78221882	BIOSAFE [®] smart220 MD				
	\Box	78221885	BIOSAFE®smart420 MD				
		78221888	BIOSAFE [®] smart500 MD				

78221894	BIOSAFE®smart1000 MD	
78221600	BIOSAFE [®] smart2400 MD	

Flyer info: Exemplary illustration. Special equipment possible.

At smart + messenger, additional costs for sms can arise. Changes and availability of the special equipment reserved.





• • • Made in Germany



BIOSAFE® VENUS

BIOSAFE® VENUS


The safe storage for your samples.



All Advantages at a Glance

The safe aluminum storage and transport vessel. Ideal for biological-medical samples.

- for up to 6,000 samples 2 ml vials in gas or liquid phase
- for up to 44,880 sequins 0,25 ml
- for up to 400 bags 25 ml






Storage capacities

BIOSAFE® VENUS aluminum series2045Liquid phase[]21,540Stainless steel canisters[piece]610Diameter canister[mm]3873Height canister[mm]120 / 280280Layers[Quantity]1 / 222 ml cryogenic vials (on vial holders)[up to]18012005 ml cryogenic vials (on vial holders)[up to]906000,25 ml sequins[up to]1.560 / 3.12016.4000,5 ml sequins[up to]720 / 1.4407.300			9	
Net capacity[I]21,540Stainless steel canisters[piece]610Diameter canister[mm]3873Height canister[mm]120 / 280280Layers[Quantity]1 / 222 ml cryogenic vials (on vial holders)[up to]18012005 ml cryogenic vials (on vial holders)[up to]906000,25 ml sequins[up to]1.560 / 3.12016.400	BIOSAFE [®] VENUS aluminum series		20	45
Stainless steel canisters[piece]610Diameter canister[mm]3873Height canister[mm]120 / 280280Layers[Quantity]1 / 222 ml cryogenic vials (on vial holders)[up to]18012005 ml cryogenic vials (on vial holders)[up to]906000,25 ml sequins[up to]1.560 / 3.12016.400	Liquid phase			
Diameter canister [mm] 38 73 Height canister [mm] 120 / 280 280 Layers [Quantity] 1 / 2 2 2 ml cryogenic vials (on vial holders) [up to] 180 1200 5 ml cryogenic vials (on vial holders) [up to] 90 600 0,25 ml sequins [up to] 1.560 / 3.120 16.400	Net capacity	[1]	21,5	40
Height canister [mm] 120 / 280 280 Layers [Quantity] 1 / 2 2 2 ml cryogenic vials (on vial holders) [up to] 180 1200 5 ml cryogenic vials (on vial holders) [up to] 90 600 0,25 ml sequins [up to] 1.560 / 3.120 16.400	Stainless steel canisters	[piece]	6	10
Layers [Quantity] 1 / 2 2 2 ml cryogenic vials (on vial holders) [up to] 180 1200 5 ml cryogenic vials (on vial holders) [up to] 90 600 0,25 ml sequins [up to] 1.560 / 3.120 16.400	Diameter canister	[mm]	38	73
2 ml cryogenic vials (on vial holders) [up to] 180 1200 5 ml cryogenic vials (on vial holders) [up to] 90 600 0,25 ml sequins [up to] 1.560 / 3.120 16.400	Height canister	[mm]	120 / 280	280
5 ml cryogenic vials (on vial holders) [up to] 90 600 0,25 ml sequins [up to] 1.560 / 3.120 16.400	Layers	[Quantity]	1/2	2
0,25 ml sequins [up to] 1.560 / 3.120 16.400	2 ml cryogenic vials (on vial holders)	[up to]	180	1200
	5 ml cryogenic vials (on vial holders)	[up to]	90	600
0,5 ml sequins [up to] 720 / 1.440 7.300	0,25 ml sequins	[up to]	1.560 / 3.120	16.400
	0,5 ml sequins	[up to]	720 / 1.440	7.300

Further vessel sizes and accessories on request.

Suitable accessories up to VENUS 45



Different colors and numbering for easier sample storage.



ful inserts Ø 65 mm



Item No.: 78216990



without lid Ø 35 mm / piece Item No.: 78212840 Ø 35 mm / Set (20 St.) Item No.: 78216987 Ø 65 mm / Stück

Item No.: 78212792



2 Lifter Item No.: 78216989



Technical Data		
BIOSAFE® VENUS		20
	-	F.0

Neck diameter	[mm]	50	120
Weight, empty	[kg]	13	24
Weight, full	[kg]	30,4	57
Diameter, outside	[mm]	388	468
Total height	[mm]	660	710
Daily evaporation	[l/day]*	0,10	0,3
Static holding time 💧	[days]*	216	144
Static holding time 👛	[days]*	_	-

*) The daily evaporation rate and the static holding time were determined at 20 degrees Celsius, 1013 hectopascal, with the vessel standing still, lid closed and without accessories. These are nominal values that may change depending on the use and age of the container.

Storage Canacities

\bigtriangledown	
	4

لمليل الماليان

Storage capacities		া			
BIOSAFE [®] VENUS aluminum series			150		180
Liquid phase / 🕮 Gas phase		6		6	
Net capacity	[]		116		172
Number of racks made of aluminum	[piece]	6	6	6	6
Number of levels	[Quantity]	8	7	10	10/9
2 ml cryogenic vials in 10x10 boxes	[up to]	4.800	4.200	6.000	5.400
Size of cryogenic boxes 10x10 (HxWxD)	[mm]	133 x 1	33 x 51	133 x ′	133 x 51
0,25 ml sequins	[up to]	33.660	22.440	44.880	33.660
0,5 ml sequins	[up to]	12.240	8.160	16.320	12.240
25 ml bags	[up to]	330	264	396	330

Further vessel sizes and accessories upon request.

Suitable accessories up to VENUS 180



Caster base

Visiotubes Item No.: 78216988

Eye protection Item No.: 78212792

Item No.: 78216985 up to VENUS 180 Item No.: 78216986

up to VENUS 45

www.cryotherm.de | www.cryotherminc.com

(for plastic canisters, available for BIOSAFE® VENUS 20).

Sample protection

A hinged lid prevents floating



All canisters are numbered, can additionally be provided with tags for individualized labels / codes. The practical hooks allow hanging onto the container neck.

Storage of sequins in Visiotubes (Ø 12 mm)



for VENUS 150/180

Cr	<u>othe</u>	erm®
	life sci	ience

450	100
150	180
215	215
42	56
136	195
586	683
962	1028
0,75	0,8
194	217
29	25



Gas phase storage

Your samples are stored in the gas phase, without contact to liquid nitrogen. Corresponding racks ensure storage in the gas phase



Liquid phase storage Your samples are stored up to 100 percent in the liquid phase at temperatures down to -196 degrees C.



Our comfortable clip closure system Every cryogenic box on each level can be easily inserted and removed.



BIOSAFE® VENUS

with telemonitoring system.



Your "Plus" for sample monitoring.



With this "PLUS" to your BIOSAFE® VENUS series you get the advantage of digital on-site monitoring of your samples. A clear, oval display offers you a monitoring system with temperature and level indicator, with RS485 interface for reports, and analog output 4-20 mAm.



BIOSAFE® VENUS COMFORT

The COMFORT version provides all the advantages of "PLUS" with the addition of an automatic refill system.



Nitrogen levels and temperature are controlled by probes and solenoids, that control the nitrogen supply of e.g. APOLLO[®].





BIOSAFE® VENUS

Order Form.

Item Numbers



For your BIOSAFE® VENUS we supply various equipment options and extensive accessories.

BIOSAFE® VENUS Vessel Options

X	Item Number		Quantit
	78216926	BIOSAFE® VENUS 20 with 6 canisters, stainless steel, 1 level	
	78216927	BIOSAFE® VENUS 20 with 6 canisters, stainless steel, 2 levels	
	78216930	BIOSAFE® VENUS 45 with 10 canisters, stainless steel, 2 levels	
	78216931	BIOSAFE® VENUS 45 PLUS with 10 canisters, stainless steel, 2 levels, level gauge and temperature display	
	78216940	BIOSAFE® VENUS 150 with 6 aluminum racks for 2 ml cryogenic vials (in liquid phase)	
	78216941	BIOSAFE® VENUS 150 PLUS with 6 aluminum racks for 2 ml cryogenic vials (in liquid phase)	
	78216942	BIOSAFE® VENUS 150 COMFORT with 6 aluminum racks for 2 ml cryogenic vials (in liquid phase)	
	78216943	BIOSAFE® VENUS 150 COMFORT with 6 aluminum racks for 2 ml cryogenic vials (in gas phase)	
	78216944	BIOSAFE® VENUS 180 with 6 aluminum racks for 2 ml cryogenic vials (in liquid phase)	
	78216945	BIOSAFE® VENUS 180 PLUS with 6 aluminum racks for 2 ml cryogenic vials (in liquid phase)	
	78216946	BIOSAFE® VENUS 180 COMFORT with 6 aluminum racks for 2 ml cryogenic vials (in liquid phase)	
	78216947	BIOSAFE® VENUS 10 COMFORT with 6 aluminum racks for 2 ml cryogenic vials (in gas phase)	

Further container sizes and accessories upon reques

BIOSAFE® VENUS Extensive accessories

X	Item Number		Quantity
	78216985	Caster base for up to BIOSAFE®VENUS 45, made of aluminum, with 3 fixation parts	
	78216986	Caster base for up to BIOSAFE®VENUS 180, made of aluminum, with 3 fixation parts	
	78212840	Transparent cup, without lid, can be ordered per item, Ø 35 mm (for BIOSAFE®VENUS 20)	
	78216987	Set of 20 transparent cups, without lid, Ø 35 mm (for BIOSAFE®VENUS 20)	
	78212792	Transparent cup, without lid, can be ordered per piece, Ø 65 m (for BIOSAFE®VENUS 45)	
	78216992	Set of 20 transparent cups, without lid, Ø 65 m (for BIOSAFE®VENUS 45)	
	78216990	Set of 5 transp. cups with 11 colorful inserts, lid, Ø 65 mm (for BIOSAFE®VENUS 45)	
	78216989	Lifters, per piece, for vessels with 2 levels	

Our extensive security equipment can be found on pages 6 and 7.



22



Our BIOSAFE® VENUS vessels are available in capacities from from 20 l to nearly 180 liters.

\$ +49 (0)2741 - 95 85 0

(f) (in) 🕑 💘 💌

BIOSAFE® MARATHON

BIOSAFE® MARATHON

For the safe shipment of your samples.



All Advantages at a Glance

The safe aluminum dry shipping container for your biological-medical samples.

• for up to 612 samples, 2 ml vials

• for up to 9,840 sequins 0,25 ml

20

20,6 1* div. div. div. 612 306

9.840

4.380

for up to 42 bags 25 ml*

● for up to 6 bags of 500 ml*



Technical Data		
BIOSAFE® MARATHON		5
Neck diameter	[mm]	50
Weight, empty, without racks	[kg]	7,5
Weight, full, without racks	[kg]	11,3
Diameter, outside	[mm]	248
Total height, without data logger	[mm]	550
Daily evaporation	[l/Tag]*	0,12
Static holding time 👛	[Tage]*	37

*) The daily evaporation rate and the static holding time were determined at 20 degrees Celsius, 1013 hectopascal, with the vessel standing still, lid closed and without accessories. These are nominal values that may change depending on the use and age of the container.

Storage capacities 🤝 🛁					
	9				
BIOSAFE® MARATHON	Aluminum	series 5	15		
Dry-Shipper					
Net capacity		6,5	15		
Number of canisters made of stainless steel	[piece]	2	2		
Diameter canister	[mm]	41	71		
Height canister	[mm]	280	280		
Levels	[quantity]	2	2		
2 ml cryogenic vials	[up to]	84	252		
5 ml cryogenic vials	[up to]	42	126		

[up to] **)Assembling bags for MARATHON 20 optional. *) Classification rack. Further vessel sizes and accessories upon request

Suitable accessories up to MARATHON 20



0,25 ml sequins

0,5 ml sequins

Caster base for MARATHON 20 Item No.: 78216985 (including 3-piece fixation set)



[up to]

1.040

400


3.280

1.320

Protective packaging-T°Tracker Item No.: 78216993 Item No.: 78216952 For your temperature Item No.: 78216994 monitoring



Cryo-Protective Shield Item No.: 78212792



¹⁾ Bag requires protective cassettes. Capacity deviations are possible depending on the manufacturer.

Cassettes with different closures sealings





www.cryotherm.de | www.cryotherminc.com

MARATHON 5


MARATHON 10


MARATHON 20

Item No.: 78216995



The BIOSAFE® MARATHON series is lockable.









BIOSAFE® MARATHON

Order form.

Item Numbers



Our BIOSAFE® MARATHONvessels are available in capacities from 5 liters to nearly 20 liters.

For BIOSAFE® VENUS we supply various equipment options and extensive accessories.

BIOSAFE[®] MARATHON Vessel options

X	Item Numbe	er Quantity
	78216949	BIOSAFE® MARATHON 5 with two canisters, stainless steel, 2 layers
	78216950	BIOSAFE® MARATHON 15 with two canisters, stainless steel, 2 layers
	78216951	BIOSAFE® MARATHON 20 with 1 rack, 5 layers, for a total of 500 samples for 2 ml vials
	78216953	BIOSAFE® MARATHON 20 with 1 rack for bags

Further vessel sizes and accessories upon request.


BIOSAFE[®] MARATHON Extensive accessories

🕅 Item Number Quant						
	78216985	Caster base for BIOSAFE® MARATHON 20, made of aluminum, with 3 fixing parts				
	78216952	Data logger for recording temperature for BIOSAFE® MARATHON series				
	78216953	Protective packaging for BIOSAFE® MARATHON 5				
	78216954	Protective packaging for BIOSAFE [®] MARATHON 15				
	78216955	Protective packaging for BIOSAFE [®] MARATHON 20				
	78216962	Rack, 5 levels, for a total of 125 samples 2ml vials (in box system 75x75x51 mm)				
	78212840	Transparent cup, without lid, can be ordered per piece, Ø 35 mm (goblet Ø 35, for BIOSAFE® VENUS 20)				
	78216987	Set of 20 transparent cups, without lid, Ø 35 mm (for BIOSAFE® VENUS 20)				
	78212792	Transparent cup, without lid, can be ordered per piece, Ø 65 m (VENUS 45)				
	78216992	Set of 20 transparent cups, without lid, Ø 65 m (for BIOSAFE®VENUS 45)				
	78216990	Set of 5 transp. cups with 11 colorful inserts, lid, Ø 65 mm (for BIOSAFE® VENUS 45)				
	78216989	Lifter, per piece, for containers with 2 layers				

• Our comprehensive safety equipment can be found on pages 6 and 7.

Sales@cryotherm.de

\$ +49 (0)2741 - 95 85 0



Unpressurized vessel series for small amounts of nitrogen.

مناطر ارارا ا



All Advantages at a Glance



Safe and lightweight aluminum storage and transport vessel for liquid nitrogen.

SMALL, COMPACT AND STURDY.

- Equipped as standard with 50 mm flange
- Available in the three most common sizes with
- a capacity of 10, 20 & 25 liters
- Extensive accessories

Technical Data

			J	
Storage Capacity ALPHA (alumi	10	20	25	
Net capacity	[]	12,2	21,5	26
Daily evaporation	[l/Day]*	0,18	0,18	0,2
Static holding time	[Days]*	67	119	130
Weight, empty	[kg]	7,5	11	13,5
Weight, full	[kg]	17,3	28,3	34,6
Total height	[mm]	584	605	669
Inner height bottom to neck	[mm]	530	535	612
Flange diameter	[mm]	50	50	50
Outer diameter	[mm]	308	388	388

*) The daily evaporation rate and the static holding time were determined at 20 degrees Celsius, 1013 hectopascal, with the vessel standing still, lid closed and without accessories. These are nominal values that may change depending on the use and age of the container.



Auxiliary handle for safe pouring of liquid nitrogen.

ALPHA Order form.

Item Numbers



For your ALPHA, different configurations are available and extensive accessories are offered.

ALPHA Vessel options


X	ltem number	
	78216920	ALPHA 10, unpressurized transport / stor
	78216921	ALPHA 20, unpressurized transport / sto
	78216922	ALPHA 25, unpressurized transport / sto

BIOSAFE® ALPHA Extensive accessories

X	Item number	
	78216985	Caster base for ALPHA 10 - 25, made of al
	78216996	Additional handle for safe pouring for ALP
	78216997	Additional handle for safe pouring for ALP
	78216998	Tipcart for ALPHA 20
	78216999	Tipcart for ALPHA 25
	78217000	Extraction attachment for low withdrawal evaporation

Our comprehensive safety equipment can be found on pages 6 and 7.

🖂 sales@cryotherm.de




Accessories available up to ALPHA 25



Auxiliary Handle For ALPHA 10 Item No.: 78216996 For ALPHA 25 Item No.: 78216997





DL3 Extraction attachment Item No.: 78217000

For ALPHA 25

Item No.: 78216999



Our BIOSAFE® ALPHA series is available with a vessel capacity of 10 to 25 liters.

Quantity

orage vessel for liquid nitrogen

orage vessel for liquid nitrogen

storage vessel for liquid nitrogen

Quantity

luminum, with 3 fixation parts

PHA 10

PHA 20 and 25

l quantities, pressure build-up by natural

C +49 (0)2741 - 95 85 0




(f) (in) 🕑 💘 💌

APOLLO[®]


Your ideal storage and transport vessel for liquid nitrogen.

APOLLO[®]

Components and flow chart.







1 Extraction valve

4 Pressure gauges

6 Level indicator

8 Phase separator

5 Safety valve

7 Filling hose

9 Castors

Essential Construction Elements and **Equipment**

مرامليلية المرامليلية

(f) (in) 🕑 💘 💌

1 Vacuum sealing and safety valve **2** Exhaust and overflow valve B Superinsolation **3** Pressure build-up valve 4 Absorbents (long-term vacuum protection) and special getters **(b** Exhaust and overflow line **1**6 Sampling line 1 Pressure build-up line 18 Chassis and protective frame Indoor and outdoor vessel Pressure regulator 20 gaseous nitrogen 2 liquid nitrogen **1** Pressure build-up regulator

APOLLO[®]

Economical. Easy-to-use. Secure.

Standard equipment

- Integrated safety device on the vessel neck
- Vacuum seal and safety device
- Smooth-running castors and transport protection frame
- Level and pressure gauge
- EK lift with exhaust and overflow valve
- Integrated pressure build-up
- Filling hose (1.5 m) with phase separator
- Transport plug **1** with safety
- Cold protection gloves 2 and safety glasses 3 with side protection



للمليلية الم

Technical Data

	APOLLO [®] 50 1,3 bar	APOLLO® 100 1,3 bar	APOLLO® 150 1,3 bar	APOLLO® 200 2,0 bar	APOLLO® 350 2,0 bar
	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel
	CE	CE	CE	CE0035	CE0035
[]	49,5	99,2	149,5	198,5	348
[bar]	1,3	1,3	1,3	2,0	2,0
[kg]	44	62	79	100	160
[kg]	85	145	204	266	450
[%/day]	2	1,2	1	0,6	0,5
[cm]	80	114	146	117	165
[cm]	50	50	50	70	70
[cm]	65	65	65	80	80
 [mm]	58,5	92	122	91	135
[mm]	125	125	125	160	160
[mm]	50	50	50	50	50
	KF 50	KF 50	KF 50	KF 50	KF 50
	[bar] [kg] [%/day] [cm] [cm] [cm] [cm] [mm] [mm]	1,3 bar Stainless steel CE [l] 49,5 [bar] 1,3 [kg] 44 [kg] [kg] 50 [cm] 65 [mm] 50 [mm]	1,3 bar 1,3 bar Stainless steel Stainless steel CE CE [1] 49,5 99,2 [bar] 1,3 1,3 [bar] 1,3 1,3 [kg] 44 62 [kg] 85 145 [%/day] 2 1,2 [cm] 80 114 [cm] 50 50 [mm] 58,5 92 [mm] 125 125 [mm] 50 50	1,3 bar 1,3 bar 1,3 bar Stainless steel Stainless steel Stainless steel Stainless steel CE CE CE [l] 49,5 99,2 149,5 [bar] 1,3 1,3 1,3 [bar] 1,3 1,3 1,3 [kg] 44 62 79 [kg] 85 145 204 [%/day] 2 1,2 1 [cm] 80 114 146 [cm] 50 50 50 [mm] 58,5 92 122 [mm] 125 125 125 [mm] 50 50 50 KF 50 KF 50 KF 50 50	1,3 bar1,3 bar2,0 barStainless steelStainless steelStainless steelStainless steelCECECECE[l]49,599,2149,5[bar]1,31,31,3[bar]1,31,32,0[kg]446279100[kg]85145204266[%/day]21,210,6[cm]80114146117[cm]50505070[cm]65656580[mm]58,59212291[mm]125125125160[mm]50505050KF 50KF 50KF 50KF 50

Upon request the APOLLO® series is also available for other types of gas.



Abgas

ТЬ

D

Entnahme

H INF 3/4"

For dimensions see table.

B

APOLLO[®]

Extensive accessories.

Options

- Pressure build-up regulator (Item No. 78202710)
- Pressure reduction regulator (Item No. 78203080)
- Double extraction liquid side (Item No. 79415416)
- Double extraction gas side (Item No. 78203081)
- Triple extraction liquid side (Item No. 79418082)
- Triple extraction gas side (Item No. 78203082)
- Safety spacer **4** (Item No. 78206124)
- Level indicator, capacitive, battery
- Forklift pockets **7**
- Set of castors, antimagnetic

Further options available upon request.



tem Numbers					
APOLLO [®] liquid nitrogen vessel	50	100	150	200	350
with analog display 5	78202700	78202701	78202702	78202748	78202749
with digital display 6	78216837	78216836	78216838	78216839	78216840
ASME Stamp 1,5 bar 🛽 🛽 🔊	78221094	78220960	78221095	78221811	78221812
Options for APOLLO [®] line					
Forklift pockets 🛛	78203083	78203083	78203083	78203975	78203900
Set of castors, antimagnetic	78200474	78200474	78200474	78203992	78203992


*) Further vessel dimensions and custom vessels available upon request.



Our APOLLO[®] series complies with the European Pressure Equipment Directive 2014/68 / EU



(f) (in) 🖻 💘 💌



www.cryotherm.de | www.cryotherminc.com







Cryotherm GmbH & Co. KG

Euteneuen 4 D-57548 Kirchen/Sieg, Germany Telephone: ++49 27 41 / 95 85-0 Fax: ++49 27 41 / 69 00 E-mail: info@cryotherm.de



Your local contact person:

Version 1119

www.cryotherm.de | www.cryotherminc.com

